



weighs only 1.75Kg making it one of the lightest syringe pumps across the globe.



**Environment** 















Q

DISPO

## **SYRU 4004**










Feather Light

Compact

Multi Porting







Uni Power Socket

Anti Bolus Function

Safety Alarm

## **SYRU 4004 Specification**

Power Supply (I/P)	100-280V AC, 50/60 Hz
Battery	12V Ni-MH Battery
Battery Backup	5hrs @ 5mL/Hr Flow rate on 50ml Syringe
Power Consumption	7 VA
Applicable Syringe	5ml, 10ml, 20ml, 30ml & 50ml Syringe from AKAS, BD, Braun, Dispovan, and Custom Programmable for any brand of Syringe.
Automatic Syringe Detector	Available
Custom Programmable Syringe	Available
Flow Rate Range	0.1 - 350 mL/Hr (5ml Syringe)
	0.1 – 500 mL/Hr (10ml Syringe)
	0.1 - 900 mL/Hr (20ml Syringe)
	0.1 - 1200 mL/Hr (30ml Syringe)
	0.1 - 2000 mL/Hr (50ml Syringe)
Infused Volume	0.1 - 9999.9 ml
Elapsed Time	00:01 to 99:59 hrs
Modes	Flow Rate only.
	Flow Rate & Target volume.
	Flow Rate & Target Time.
	Target volume and target time
	High Pressure mode
Bolus	Purge Bolus & Auto Bolus With separate key access for quick infusion
Bolus Flow Rate	0.1 - 350 mL/Hr (5ml Syringe)
	0.1 – 500 mL/Hr (10ml Syringe)
	0.1 - 900 mL/Hr (20ml Syringe)
	0.1 - 1200 mL/Hr (30ml Syringe)
	0.1 - 2000 mL/Hr (50ml Syringe)
Anti Bolus Function	Post occlusion – Bolus reduction to avoid accidental bolus infusion
Bolus Volume	0.1 mL to maximum volume of the syringe loaded.
KVO	0.1mL/Hr to 20mL/Hr in 0.1 ml increment.
Occlusion Pressure Range	100 – 1000 mmHg in increments of 50mmHg.

Alarms	AL 1	Syringe holder disturbed during infusion Motor unable to drive	
	AL 3	Priming not done Occlusion sensor Failed Plunger not engaged	
	AL 10	Drive Disengaged Syringe Load error Target reached End of	
	AL 12	Syringe Occlusion Pre Occlusion Warning Battery discharged	
	AL 13	Target, End, Near Empty, Battery Low, Operation on mains / battery, Bolus, KVO, Occlusion level, Prime, Near to Target	
	AL 14 tAr.	LED Display	
	End	Mechanical Accuracy ±1% Including Syringe ± 2%	
	occL	Adjustable 5 levels Possible to change the infusion rate without	
	P.Ocl OFF	stopping the infusion Save last infusion data even when power on	
	OFF	& Off	
		Dock 4 pumps & mountable on poles with single power supply-chord $$	
		170 X 107 X 135 ( L x B x H) 2.2 ltrs 1.75 kg	
Audio/Visual Indication		ISO 13485, EUROPEAN CE, CDSCO registered Medical Device Manufacturer BF type, Class II protection rating, in compliance	
Display		with IEC/EN 60601-1, Class IIb in compliance with Council Directive 93/42/EEC Temperature: -10 to 50 'C Humidity: 15 to 90%	
Accuracy		Pressure: 0 to 106 KPa EN 60601-1, EN 60601-1-2, EN 60601-1-8, EN 60601-2-24	
Alarm LED Brightness		EN 00001-1, EN 00001-1-2, EN 00001-1-8, EN 00001-2-24	
Inline changing of flowrat	e		
Last infusion data			
Docking Station			
Dimension			
Space Occupied			
Weight			
Certification			
Classification			
Operating & Storage Cond	ditions		
Safety Testing			